

Glass DC Nebulizers

Conikal
1mL/min uptake A41-1-UC1S

SeaSpray
400µL/min uptake A41-07-USS04

Slurry
4mL/min uptake A41-1-US6

MicroMist
600µL/min uptake A41-07-UM06S
200µL/min uptake A41-07-UM02
100µL/min uptake A41-07-UM01
50µL/min uptake A41-07-UM005



DC Nebulizers Compatible with HF

OpalMist
1mL/min uptake A41-1-PFA1S
400µL/min uptake A41-07-PFA04
200µL/min uptake A41-07-PFA02
100µL/min uptake A41-07-PFA01
50µL/min uptake A41-07-PFA005



DuraMist
1mL/min uptake A41-1-DM1S
400µL/min uptake A41-07-DM04



VeeSpray
0.6 - 3mL/min uptake A41-1-CV6



Spray Chambers

Glass
50mL Cyclonic 21-809-3994
With B14 gravity drain



Quartz
50mL Cyclonic 21-809-3992
With B14 gravity drain



Drain Device 21-808-8214



Extention Pipe
Glass Extention w/ right angle 31-808-3920

Spray Chambers Compatible with HF

PFA
50mL Cyclonic 21-809-3991



HF Drain Device 21-809-2540




Extention Pipe
PFA Extention w/ right angle 21-809-4081

Torches

Quartz Torches

| | |
|---|-------------|
| Quartz Mini Torch (Suitable for ICP-MS 2030) | 30-807-0567 |
| Torch with 1.2mm Injector (Suitable for ICP-MS 8500) | 30-808-9986 |
| Torch with 1.5mm Injector (Suitable for ICP-MS 8500) | 30-808-1090 |



RF Coils

| | |
|--------|--------------|
| Copper | 70-900-6100C |
| Silver | 70-900-6100S |
| Gold | 70-900-6100G |

RF Coil Installation Kit

| | |
|--------------------------|--------------|
| RF Coil Installation Kit | 70-900-6100T |
|--------------------------|--------------|

Cones

Sampler Cones

| | |
|------------------|--------------|
| Nickel Sampler | SZ7001-Ni |
| Copper Sampler | SZ7001-Cu |
| Platinum Sampler | SZ7003-Pt/Cu |

Skimmer Cones

| | |
|------------------|--------------|
| Copper Skimmer | SZ7002-Cu |
| Nickel Skimmer | SZ7002-Ni |
| Platinum Skimmer | SZ7004-Pt/Cu |

Common Accessories

Elegra Argon Humidifier

| | |
|-------------------------|-------------|
| Elegra Argon Humidifier | 70-803-1312 |
|-------------------------|-------------|



Gaurdian In-Line Filters

| | |
|---|-------------|
| Guardian in-line non-return gas filter | 70-803-1362 |
| Guardian in-line sample filter | 70-803-1108 |



Trident Internal Standard Kit

| | |
|---------------------|-------------|
| Standard Kit | 60-808-1185 |
| Standard Kit for HF | 60-808-1150 |



Autosampler Probes

(See page 12 for complete list of probes)

| | |
|--|-------------|
| Autosampler Probe for Shimadzu ASC-6100F/9800 | 70-803-1200 |
|--|-------------|