

ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

ISSUE 9

Product List for Agilent (Varian) ICP-OES: 700-ES Series Axial, Vista Axial



GLASS EXPANSION
Quality By Design

INTERNATIONAL



15 Batman Street
West Melbourne
Vic 3003, Australia


Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com



AMERICAS

4 Barlows Landing Road, Unit 2A
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097
Telephone: 508 563 1800
Facsimile: 508 563 1802
Email: geusa@geicp.com

Glass Nebulizers	
Conikal 2mL/min uptake ARG-07-UC2	
SeaSpray 2mL/min uptake ARG-07-USS2	
Slurry 4mL/min uptake ARG-07-US6	
MicroMist 600µL/min uptake ARG-07-UM06 200µL/min uptake ARG-07-UM02 100µL/min uptake ARG-07-UM01 50µL/min uptake ARG-07-UM005	
Nebulizers suitable for use with HF	
OpalMist 2mL/min uptake ARG-07-PFA2 600µL/min uptake ARG-07-PFA06 200µL/min uptake ARG-07-PFA02 100µL/min uptake ARG-07-PFA01 50µL/min uptake ARG-07-PFA005	
DuraMist 1mL/min uptake ARG-07-DM1 400µL/min uptake ARG-07-DM04	
VeeSpray - Ceramic Standard ARG-07-CV6	

IsoMist	
IsoMist Programmable Temperature Spray Chamber for 700-ES or Vista Axial	
With glass spray chamber KT-1014	
With quartz spray chamber KT-1014Q	
With PFA spray chamber KT-1014P	
For all glass system KT-1054	

Spray Chambers		
Tracey 50mL Cyclonic 20-809-0600HE		
Twister 50mL Cyclonic 20-809-0285HE		
Cinnabar 20mL Cyclonic 20-809-0183HE		
Twinnabar 20mL Cyclonic 20-809-0448HE		
Tracey TFE 50mL Cyclonic 20-809-9188 Suitable for use with HF.		
Mounting Bracket - HF Spray Chamber 21-100-0364		
PuraFlex Transfer Tube 21-809-2831 Connects spray chamber to torch. Requires adjustable mount 21-100-0215, or for HF, 21-100-0364.		

Agilent (Varian) ICP-OES: 700-ES Series Axial, Vista Axial

Torches

D-Torch

[30-808-2919](#)

Comprises:

Base and Inner Tube [31-808-2918](#)

Quartz Outer Tube [31-808-2685](#)

Retaining Ring for Outer Tube [31-808-2818](#)

Retaining Ring for Injector [31-808-2921](#)

Options:

Tapered Quartz Injector 2.4mm [31-808-3043](#)

Tapered Alumina Injector 0.8mm [31-808-3364](#)

Alumina Injector 2.4mm [31-808-3133](#)

Alumina Injector 1.8mm [31-808-3325](#)

Sapphire Injector 2.0mm [31-808-3355](#)

Ceramic Outer Tube [31-808-3027](#)

PuraFlex Transfer Tube with bend [21-809-3707](#)

Fixed Torches


Quartz Torch with 90° Bend and 2.3mm Injector [30-808-0478](#)

Quartz Torch (High Solids, short) with 90° Bend and 2.3mm Injector [30-808-9983](#)

Quartz Torch (High Solids, Full-Length) with 90° Bend and 2.3mm Injector [30-807-0528](#)

Quartz Torch with 90° Bend and 1.8mm Injector [30-808-3315](#)

Quartz Torch with 90° Bend and 0.8mm Injector [30-808-1346](#)



Torches (continued)

Semi Demountable Torch

[30-808-0379](#)

Comprises:

Quartz Torch Body [31-808-0353](#)

Tapered Quartz Injector 2.4mm [31-808-0400](#)

Torch Adaptor [31-808-0863](#)

GazFit Connectors (PKT 4) - for 6mm OD side arms [GAZ-06](#)

GazFit Connectors (PKT 2) - for 8mm torch adaptor [GAZ-08](#)

Semi Demountable Torch Options

Quartz Torch Body - Short (for high solids) [31-808-0177](#)

Tapered Quartz Injector 1.0mm [31-808-0544](#)

Tapered Quartz Injector 1.5mm [31-808-0543](#)

Tapered Quartz Injector 1.8mm [31-808-0545](#)

Tapered Quartz Injector 2.4mm [31-808-0400](#)




Fixed Torches for All Glass System

Quartz Torch with 2.3mm Injector and Right Angle Ball Joint [30-808-0540](#)

Quartz Torch with 1.8mm Capillary Injector - (for oil analysis) and Right Angle Ball Joint [30-808-0887](#)

Quartz Torch (High Solids) with 2.3mm Injector and Right Angle Ball Joint [30-808-1255](#)

Quartz Torch (organics) with 0.8mm Injector and Right Angle Ball Joint [30-808-8216](#)



Semi Demountable Torch with Alumina Injector

[30-808-0383](#)

Comprises:

Quartz Torch Body [31-808-0353](#)

Tapered Alumina Injector 2.4mm [31-808-0401](#)

Torch Adaptor [31-808-0863](#)

GazFit Connectors (PKT 4) - for 6mm OD side arms [GAZ-06](#)

GazFit Connectors (PKT 2) - for 8mm torch adaptor [GAZ-08](#)



HF Resistant Injectors for Semi Demountable Torch









Tapered Alumina Injector 1.8mm [31-808-0814](#)

Tapered Alumina Injector 2.4mm [31-808-0401](#)

Sapphire Injector 2.4mm [31-808-1426](#)



RF Coils

<p>RF Coil Gold Vista 70-900-1000G</p>	
<p>RF Coil Silver Vista 70-900-1000S</p>	
<p>RF Coil Silver, Teflon Coated Vista 70-900-1000ST</p>	
<p>Installation kit Vista 70-900-1000T</p>	
<p>RF Coil Gold 700ES 70-900-1004G</p>	
<p>RF Coil Silver 700ES 70-900-1004S</p>	
<p>RF Coil Silver, Teflon Coated 700ES 70-900-1004ST</p>	
<p>Installation kit 700ES 70-900-1004T</p>	

Kits

<p>Axial All Glass System KT-8045A</p>	
<p>Kit Comprises:</p>	
<p>Twister Spray Chamber 20-809-9199HE</p>	
<p>Quartz Torch RA/FS13 & 2.3mm Injector 30-808-0540</p>	
<p>Glass Transfer Tube-206, + JC-13 Clip 31-809-0274</p>	
<p>Adjustable Mount 21-100-0215</p>	