

ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

Issue 9

Product List for Agilent (Varian) ICP-OES: 700-ES Series Radial, Vista Radial



GLASS EXPANSION
Quality By Design

INTERNATIONAL

15 Batman Street
West Melbourne
Vic 3003, Australia

Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com

AMERICAS

4 Barlows Landing Road, Unit 2A
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097
Telephone: 508 563 1800
Facsimile: 508 563 1802
Email: geusa@geicp.com

Agilent (Varian) ICP-OES: 700-ES Series Radial, Vista Radial

Glass Nebulizers

Conikal
2mL/min uptake [ARG-07-UC2](#)

SeaSpray
2mL/min uptake [ARG-07-USS2](#)

Slurry
4mL/min uptake [ARG-07-US6](#)

MicroMist
600µL/min uptake [ARG-07-UM06](#)
200µL/min uptake [ARG-07-UM02](#)
100µL/min uptake [ARG-07-UM01](#)
50µL/min uptake [ARG-07-UM005](#)



Nebulizers suitable for use with HF

OpalMist
2mL/min uptake [ARG-07-PFA2](#)
600µL/min uptake [ARG-07-PFA06](#)
200µL/min uptake [ARG-07-PFA02](#)
100µL/min uptake [ARG-07-PFA01](#)
50µL/min uptake [ARG-07-PFA005](#)



DuraMist
1mL/min uptake [ARG-07-DM1](#)
400µL/min uptake [ARG-07-DM04](#)



VeeSpray - Ceramic
Standard [ARG-07-CV6](#)



IsoMist

IsoMist Programmable Temperature Spray Chamber for 700-ES or Vista Radial

With glass spray chamber [KT-1022](#)
With quartz spray chamber [KT-1022Q](#)
With PFA spray chamber [KT-1022P](#)
For all Glass System [KT-1055](#)




Spray Chambers

Tracey
50mL Cyclonic [20-809-0600HE](#)



Twister
50mL Cyclonic [20-809-0285HE](#)



Cinnabar
20mL Cyclonic [20-809-0183HE](#)



Twinnabar
20mL Cyclonic [20-809-0448HE](#)



Tracey TFE
50mL Cyclonic [20-809-9188](#)
Suitable for use with HF.



Mounting Bracket - HF Spray Chamber
[21-809-2555](#)



PuraFlex Transfer Tube
[21-809-2735](#)
Connects spray chamber to torch.
Requires adjustable mount 21-809-2555.



Torches

D-Torch Demountable Torch

D-Torch with Quartz Outer Tube [30-808-3052](#)

Comprises:

Base and Inner Tube [31-808-3146](#)

Quartz Outer Tube [31-808-3062](#)

Retaining Ring for Outer Tube [31-808-2818](#)

Retaining Ring for Injector [31-808-2921](#)

Options:

Tapered Quartz Injector 1.5mm [31-808-3058](#)

Tapered Quartz Injector 0.8mm [31-808-3195](#)

Ceramic Outer Tube [31-808-3092](#)

Tapered Quartz Injector 1.5mm with Ball Joint [31-808-3150](#)

Tapered Quartz Injector 2.4mm with Ball Joint [31-808-2925](#)

Alumina Injector 1.8mm [31-808-3160](#)

Alumina Injector 1.5mm [31-808-3230](#)

Sapphire Injector 2.0mm [31-808-3312](#)



Semi Demountable Torch

[30-808-0542](#)

Comprises:

Torch Body, Quartz [31-808-0541](#)

Torch Adaptor 8/44 [31-808-2638](#)

Tapered Quartz Injector 1.5mm [31-808-0543](#)

GazFit Connectors (PKT 4) for 6mm OD side arms [GAZ-06](#)

Quartz Injectors

Tapered Quartz Injector 1.0mm [31-808-0544](#)

Tapered Quartz Injector 1.5mm [31-808-0543](#)

Tapered Quartz Injector 1.8mm [31-808-0545](#)

Tapered Quartz Injector 2.4mm [31-808-0400](#)



Fixed Torches

Quartz Torch, Low Flow, 1.5mm Injector [30-800-0015](#)

Quartz Torch, High Solids, 1.8mm Injector [30-808-1182](#)

Quartz Torch, 0.8mm Injector [30-807-0539](#)



Other Products

Quartz Torch Bonnet [31-800-0003](#)



Torches (continued)

Semi Demountable Torch HF

[30-808-0942](#)

Comprises:

Torch Body, Quartz [31-808-0541](#)

Torch Adaptor 8/44 [31-808-2638](#)

Tapered Alumina Injector 1.8mm [31-808-0938](#)

GazFit Connectors (PKT 4) for 6mm OD side arms [GAZ-06](#)



RF Coils

RF Coil

Gold Vista [70-900-1001G](#)



RF Coil

Silver Vista [70-900-1001S](#)



RF Coil

Silver, Teflon Coated Vista [70-900-1001ST](#)



Installation kit

Vista [70-900-1001T](#)



RF Coil

Gold 700ES [70-900-1005G](#)



RF Coil

Silver 700ES [70-900-1005S](#)



RF Coil

Silver, Teflon Coated 700ES [70-900-1005ST](#)



Installation kit

700ES [70-900-1005T](#)



Kits

Radial All Glass System Kit

[KT-8046A](#)

Kit Comprises:

Twister Spray Chamber [20-809-9199HE](#)

Radial Torch with Ball Joint and 1.5mm Injector [30-808-0875](#)

Glass Transfer Tube-160, + JC-13 Clip [31-809-0549](#)

Adjustable Mount [21-100-0200](#)

All Glass System Options

Jacketed Spray Chamber Adjustable Mount [21-100-0217](#)

Quartz Torch with Ball Joint and 1.8mm Injector High Solids [30-807-0536](#)

Quartz Torch with Ball Joint and 2.4mm Injector High Solids [30-807-0569](#)

