

ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

Issue 9

Product List for PerkinElmer ICP-OES: Optima 4300V, 5300V, 7300V



INTERNATIONAL

15 Batman Street
West Melbourne
Vic 3003, Australia

Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com

AMERICAS

4 Barlows Landing Road, Unit 2A
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097
Telephone: 508 563 1800
Facsimile: 508 563 1802
Email: geusa@geicp.com

PerkinElmer ICP-OES: Optima 4300V, 5300V, 7300V

Glass Nebulizers

Conikal	
2mL/min uptake	ARG-07-UC2
SeaSpray	
2mL/min uptake	ARG-07-USS2
Slurry	
4mL/min uptake	ARG-07-US6
MicroMist	
600µL/min uptake	ARG-07-UM06
200µL/min uptake	ARG-07-UM02
100µL/min uptake	ARG-07-UM01
50µL/min uptake	ARG-07-UM005



Nebulizers suitable for use with HF

OpalMist	
2mL/min uptake	ARG-07-PFA2
600µL/min uptake	ARG-07-PFA06
200µL/min uptake	ARG-07-PFA02
100µL/min uptake	ARG-07-PFA01
50µL/min uptake	ARG-07-PFA005
DuraMist	
1mL/min uptake	ARG-07-DM1
400µL/min uptake	ARG-07-DM04
VeeSpray - Ceramic	
Standard	ARG-07-CV6



IsoMist


IsoMist Programmable Temperature Spray Chamber for PerkinElmer Optima 4300/5300/7300V

With glass spray chamber	KT-1011
With quartz spray chamber	KT-1011Q
With PFA spray chamber	KT-1011P



Spray Chambers

Twister	
50mL Cyclonic	20-809-0495HE
Tracey TFE	
50mL Cyclonic	20-809-2507
Suitable for use with HF.	



Torches


D-Torch Demountable Torch	
D-Torch for 4300/5300/7300V	30-808-2832
Comprises:	
Base and Inner Tube	31-808-2809
Quartz Outer Tube	31-808-2821
Retaining Ring	31-808-2818



Options:	
Ceramic Outer Tube	31-808-2815



Quartz Torch	
Quartz Torch - three slot	30-807-0064



Quartz Injectors

Tapered Quartz Injector 1.0mm	31-808-0119
Tapered Quartz Injector 1.5mm	31-808-0124
Tapered Quartz Injector 1.8mm	31-807-0008
Tapered Quartz Injector 2.0mm	31-807-0050

These injectors are suitable for both the fixed torch and D-Torch.



HF Resistant Injectors

Tapered Alumina Injector 1.5mm	31-807-0051
Tapered Alumina Injector 1.8mm	31-807-0052
Tapered Alumina Injector 2.0mm	31-808-1352
Capillary Alumina Injector 2.0mm	31-808-2625

These injectors are suitable for both the fixed torch and D-torch.




Torches (continued)

Torch Bonnet

[31-800-0018](#)Purge Extension Window
for 4300/5300/7300V[31-808-1345](#)

Torch Adaptor with Auxiliary Port

[31-808-2503](#)

RF Coils

RF Coil

Gold

[70-900-2003G](#)

RF Coil

Silver

[70-900-2003S](#)

Installation Kit

[70-900-2003T](#)