

ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

Issue 9

Product List for Spectro (Ametek) ICP-OES: EOP (Modula, CIROS, Genesis, ARCOS)



GLASS EXPANSION
Quality By Design

INTERNATIONAL

15 Batman Street
West Melbourne
Vic 3003, Australia

Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com

AMERICAS

4 Barlows Landing Road, Unit 2A
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097
Telephone: 508 563 1800
Facsimile: 508 563 1802
Email: geusa@geicp.com

Spectro (Ametek) ICP-OES: EOP (Modula, CIROS, Genesis, ARCOS)

Glass Nebulizers

Conikal
2mL/min uptake [ARG-07-UC2](#)

SeaSpray
2mL/min uptake [ARG-07-USS2](#)

Slurry
4mL/min uptake [ARG-07-US6](#)

MicroMist
600µL/min uptake [ARG-07-UM06](#)
200µL/min uptake [ARG-07-UM02](#)
100µL/min uptake [ARG-07-UM01](#)
50µL/min uptake [ARG-07-UM005](#)



Nebulizers suitable for use with HF

OpalMist
2mL/min uptake [ARG-07-PFA2](#)
600µL/min uptake [ARG-07-PFA06](#)
200µL/min uptake [ARG-07-PFA02](#)
100µL/min uptake [ARG-07-PFA01](#)
50µL/min uptake [ARG-07-PFA005](#)



DuraMist
1mL/min uptake [ARG-07-DM1](#)
400µL/min uptake [ARG-07-DM04](#)



VeeSpray - Ceramic
Standard [ARG-07-CV6](#)



IsoMist

IsoMist Programmable Temperature Spray Chamber

With glass spray chamber [KT-1027](#)
With quartz spray chamber [KT-1027Q](#)
With PFA spray chamber [KT-1027P](#)



Spray Chambers

Tracey
50mL Cyclonic [20-809-0439HE](#)



Spray Chambers (continued)

Twister
50mL Cyclonic [20-809-0243HE](#)



Tracey TFE
50mL Cyclonic [20-809-2506](#)
Suitable for use with HF.
Requires Transfer Tube 21-809-2582.




PuraFlex Transfer Tube
[21-809-2582](#)



Torches

D-Torch with Quartz Outer Tube
[30-808-3371](#)

Comprises:
Base and Inner Tube [31-808-3372](#)
Quartz Outer Tube [31-808-3142](#)
Injector Adaptor [31-808-3374](#)
Injector Ferrule [31-808-2962](#)
Retaining Ring [31-808-2818](#)



Options:
Ceramic Outer Tube (requires fibre optic below) [31-808-3250](#)



Fibre Optic [31-808-3380](#)



Tapered Quartz Injector 1.8mm [31-808-2974](#)



Tapered Alumina Injector 1.2mm [31-808-3238](#)



Tapered Alumina Injector 1.8mm [31-808-2963](#)



Tapered Alumina Injector 2.4mm [31-808-3388](#)



Sapphire Injector 2.4mm [31-808-3414](#)




Torches (continued)

Quartz Torches

EOP Quartz Torch
with 2.5mm Injector [30-808-0317](#)



EOP Fixed Sheath Gas Torch
with 3.0mm Injector [30-808-1165](#)



EOP Quartz Torch
with 1.8mm Injector [30-808-8005](#)

**Semi Demountable Torch (HF)**

Semi Demountable Torch
for Spectro EOP [30-808-1467](#)

Comprises:

SDT Quartz Torch Body [31-808-1458](#)

Tapered Alumina Injector 2.4mm [31-808-1195](#)

SDT Adaptor [31-808-1459](#)

