

### Glass DC Nebulizers

**SeaSpray**

2mL/min uptake      A13-07-USS2

**MicroMist**

400µL/min uptake      A13-05-UM04

**Conikal**

1mL/min uptake      A13-05-UC1



### DC Nebulizers for HF

**DuraMist**

1mL/min uptake      A13-05-DM1

**OpalMist**

1mL/min uptake      A13-05-PFA1



### IsoMist XR

**Programmable Temperature Cyclonic Spray Chamber**

IsoMist XR Kit      KT-1135-XR

**Comprises:**

Glass Twister Spray Chamber      21-809-3185

Assembly (LOL)      20-809-5175

PuraFlex Transfer Tube (PQ9200)      21-809-5273

PuraFlex Transfer Tube      21-809-2962

Metal Ball Joint Clip      JC-13

UniFit Right Angle Connector (PKT 3)      GAZ-04-90



### BC Spray Chambers

**Tracey BC**

30mL Cyclonic      20-809-5318

\* Direct replacement for 20-808-8882HE

**NEW\***



**Tracey BC PEEK**

30mL Cyclonic      20-809-5300

\* Direct replacement for 20-809-2506

\* No internal resurfacing required

**NEW\***



### Spray Chambers

**HydraMist Simultaneous Cold Vapor/Pneumatic Nebulization Spray Chamber**

Hydride Kit      KT-1157



**Twister**

50mL Cyclonic      20-809-0222HE



**Tracey**

50mL Cyclonic      20-808-8882HE



**Cinnabar**

20mL Cyclonic      20-809-0164HE



### Spray Chambers for HF

**Tracey TFE Spray Chamber**

50mL Cyclonic      20-809-2506



## Torches

### Fully Demountable Torch

Torch Body 31-808-8896



### Fully Demountable Torch Options:

Quartz Outer Tube 31-808-8020



Syalon Outer Tube 31-808-0023

Quartz Inner Tube 31-808-8145



Alumina Inner Tube 31-808-0017

Ball Joint Adaptor 31-807-0018

Torch Bonnet - Large 31-808-3607

Radial Purge Window 31-808-8098

O-rings for Glass Body (PKT 10) 70-4X2.5V75

Viton O-rings for FDT/SDT Ball Joints (PKT 10) 70-V-011

PTFE Coated Viton O-rings for Cyclonics/Nebulizer Seal (PKT 10) VT-009

EzyGlide Cloth 70-EZG-CLOTH



### Tapered Injectors

Quartz Injector 1.0mm 31-808-0009



Quartz Injector 2.0mm 31-808-0008

Alumina Injector 2.0mm 31-808-0027



## Common Accessories

### Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1266



### Guardian In-Line Filters

Guardian in-line sample filter 70-803-1108

Eluo Adaptor for Inline Filter 70-803-1160




## Common Accessories

(continued)

### Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



### Magnifier Inspection Tool 70-803-1923



## Autosampler Probes (See page 22-23 for complete list of probes)

### For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793



### For Cetac ASX-7400, 7600 Series

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-2940



### For SPS 3/SPS 4

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-1957

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-0908

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0853

