

## Glass DC Nebulizers

### SeaSpray

400µL/min uptake      A23-1-USS04

### MicroMist

400µL/min uptake      A23-1-UM04

### Conikal

1mL/min uptake      A23-1-UC1



## Quartz DC Nebulizers

### SeaSpray

1mL/min uptake      A23-1-QSS1



## DC Nebulizers for HF

### DuraMist

400µL/min uptake      A23-1-DM04



### OpalMist

400µL/min uptake      A23-1-PFA04



## Spray Chambers

### Twister

50mL Cyclonic      20-809-4417HE

50mL Cyclonic      20-809-4560

With Matrix Port



### Tracey

50mL Cyclonic      20-809-4419HE



## Spray Chambers


### Cinnabar

20mL Cyclonic      20-809-4480HE



### High Efficiency Sample Introduction System

HE-SIS for PerkinElmer NexION      KT-1184



## Torches

### D-Torch

30-808-3927

#### Comprises:

- Base and Inner Tube      31-808-3928
- Quartz Outer Tube      31-808-3942
- Retaining Ring      31-808-2818
- Injector Adaptor      31-808-3980
- Injector Ferrule      31-808-2962
- Fastener Assembly      31-808-4130
- Viton O-rings (PKT 10)      70-V-011
- PTFE Coated Viton O-rings (PKT 10)      70-VT-012



### Quartz Injectors

Tapered 1.0mm      31-808-3423

Tapered 1.8mm      31-808-4005

Tapered 2.0mm      31-808-4006

Tapered 2.5mm      31-808-3302



### HF Resistant Injectors

Tapered Alumina 2.4mm      31-808-3205

Sapphire 2.4mm      31-808-3378



## Cones

### Sampler Cones

Nickel Sampler PE3011-Ni



Platinum Sampler PE3013-Pt

Platinum Sampler (18mm) PE3013A-Pt

Platinum Sampler (12mm) PE3013C-Pt

Platinum Sampler PE3013-Pt

Boron Free PE3013-Pt-BF



### Skimmer Cones

Nickel Skimmer PE4212-Ni

Platinum Skimmer PE4214-Pt

Nickel Hyper Skimmer PE4015-Ni

All Nickel Hyper Skimmer PE4016-Ni



### O-ring for NexION Hyper Skimmer Cone (pack of 5)

PE5125



### ConeGuard Thread Protectors

Sampler Cone Guard 70-803-1024  
For NexION



### Gasket for Sampler Cone (pack of 5)

PE5111



## Common Accessories

### Elegra Argon Humidifier

Elegra Dual Argon Humidifier 70-803-1473



### Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



### Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



### Magnifier Inspection Tool 70-803-1923



## Autosampler Probes (See page 22-23 for complete list of probes)

### For PerkinElmer S20 Series

Guardian Probe 0.75mm  
Probe Connecting Line (Red) 70-803-2097



### For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm  
Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler  
Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre  
Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre  
Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre  
Probe 1.0mm ID 70-803-0793

