

Glass DC Nebulizers

| | |
|-----------------|-------------|
| Conikal | |
| 1mL/min uptake | A41-1-UC1S |
| SeaSpray | |
| 2mL/min uptake | A41-1-USS2S |
| Slurry | |
| 4mL/min uptake | A41-1-US6S |



DC Nebulizers for HF

| | |
|-----------------|-------------|
| OpalMist | |
| 1mL/min uptake | A41-1-PFA1S |
| DuraMist | |
| 1mL/min uptake | A41-1-DM1S |



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

| | |
|-----------------------------|-------------|
| IsoMist XR Kit | KT-1067-XR |
| Comprises: | |
| Glass Twister Spray Chamber | |
| With B14 drain | 21-809-3270 |
| Assembly B14 (LOL) | 20-809-5187 |
| Adaptor | 21-809-3998 |
| Adjustable Stand | 70-803-2798 |
| Metal Ball Joint Clip | JC-13 |
| Drain Device | 21-808-8214 |



Programmable Temperature Cyclonic Spray Chamber

| | |
|-----------------------------|-------------|
| IsoMist XR Kit | KT-2749-XR |
| Requires peristaltic pump | |
| Comprises: | |
| Glass Twister Spray Chamber | 21-809-3186 |
| Assembly (ROL) | 20-809-5176 |
| Adaptor | 21-809-3998 |
| Adjustable Stand | 21-809-4911 |
| Metal Ball Joint Clip | JC-13 |
| UniFit Drain (PKT 3) | UFT-16-75 |



Spray Chambers

| | |
|------------------------|---------------|
| Tracey | |
| 50mL Cyclonic | 20-808-8881HE |
| With B14 gravity drain | |



| | |
|------------------------|---------------|
| Twister | |
| 50mL Cyclonic | 20-809-9119HE |
| With B14 gravity drain | |



| | |
|----------------|---------------|
| Twister | |
| 50mL Cyclonic | 20-809-0222HE |



| | |
|-----------------|---------------|
| Cinnabar | |
| 20mL Cyclonic | 20-809-0164HE |



HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

| | |
|-------------|---------|
| Hydride Kit | KT-1157 |
|-------------|---------|



Spray Chambers

(continued)

Drain Device 21-808-8214



Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch

Glass Transfer Tube 86mm 21-809-2846



Spray Chambers for HF

Tracey TFE

50mL Cyclonic 20-809-2508

With B14 gravity drain
Suitable for use with HF



HF Drain Device 21-809-2540



Mounting bracket for HF Spray Chamber

21-809-2839



Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch

PFA Transfer Tube 82.5mm 31-808-3071



Torches

Semi Demountable Torch (HF)

30-808-4335

Suitable for ICPS-7500/8100 and ICPE-9000/9800

Comprises:

Torch Body 31-808-4341

PEEK Injector Adaptor 31-808-4478

Tapered Alumina Injector 1.5mm 31-808-3202

Torch Options:

Tapered Quartz Injector 1.0mm 31-808-3423

Tapered Quartz Injector 1.8mm 31-808-4005

Tapered Quartz Injector 2.0mm 31-808-4006

Tapered Quartz Injector 2.5mm 31-808-3302

Tapered Alumina Injector 2.4mm 31-808-3205

Sapphire Injector 2.4mm 31-808-3378

Viton O-rings for FDT/SDT ball joints (PKT 10) 70-V-011

GazFit Connectors (PKT 4) for 4mm OD side arms GAZ-04

Quartz Torches


Quartz Torch Standard 1.5mm Injector (Suitable for 7500/8100) 30-808-1179

Quartz Torch for Organics, 1.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0519

Quartz Torch for High TDS, 1.8mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0518

Quartz Torch for High TDS, 2.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0564

Quartz Mini Torch 1.5mm Injector (Suitable for ICPE-9000/9800) 30-807-0534



RF Coils

| | |
|--------|--------------|
| Copper | 70-900-6100C |
| Silver | 70-900-6100S |
| Gold | 70-900-6100G |



RF Coil Installation Kit 70-900-6100T



Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1312




Guardian In-Line Filters

| | |
|---|-------------|
| Guardian DC in-line non-return gas filter | 70-803-1942 |
| Guardian in-line sample filter | 70-803-1108 |
| Eluo Adaptor for Inline Filter | 70-803-1160 |



Trident Internal Standard Kits

| | |
|-------------------------|-------------|
| Trident CT Standard Kit | 60-703-1179 |
| Standard Kit for HF | 60-808-1150 |



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 22-23 for complete list of probes)

For Cetac ASX-200/500/800 Series

| | |
|--------------------------------------|--------------|
| Guardian Probe 0.75mm | |
| Probe Connecting Line (Red) | 70-803-1803 |
| Polyimide sheathed Autosampler Probe | |
| 0.5mm ID | 60-808-1186L |
| PTFE Encapsulated Carbon Fibre Probe | |
| 0.5mm ID | 70-803-0784 |
| PTFE Encapsulated Carbon Fibre Probe | |
| 0.75mm ID | 70-803-1443 |
| PTFE Encapsulated Carbon Fibre Probe | |
| 1.0mm ID | 70-803-0793 |



For Cetac ASX-7400, 7600 Series

| | |
|-----------------------------|-------------|
| Guardian Probe 0.75mm | |
| Probe Connecting Line (Red) | 70-803-2940 |



For Shimadzu AS-10 & AS-20

| | |
|----------------------------------|-------------|
| Guardian Probe 0.75mm | |
| Probe Connecting Line (Red) | 70-803-2872 |
| PTFE Sheathed Carbon Fibre Probe | |
| 0.75mm ID with UniFit | 70-803-1477 |

