

## Glass DC Nebulizers

### SeaSpray

2mL/min uptake      A21-1-USS2

### MicroMist

400µL/min uptake      A21-1-UM04

### Slurry

4mL/min uptake      A21-1-US6

### Conikal

2mL/min uptake      A21-1-UC2



## Spray Chambers

### Tracey

50mL Cyclonic      20-809-2665HE

### Twister

50mL Cyclonic      20-809-2951HE

### Scott Spray Chamber

20-809-4053



## DC Nebulizers for HF

### DuraMist

1mL/min uptake      A21-1-DM1



### HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

Hydride Kit      KT-1160



## IsoMist XR

### Programmable Temperature Cyclonic Spray Chamber

IsoMist XR Kit      KT-1141-XR

#### Comprises:

Glass Twister Spray Chamber      21-809-3185

Assembly (LOL)      20-809-5175

Adaptor 12/8 for SpectroBlue      31-808-3489

Glass Ball Cup Adaptor S19      21-809-4101

Tygon Tubing 1/4 inch x 1m      1921-1M

Metal Ball Joint Clip      JC-19

GazFit Connectors for 8mm  
OD arm (PKT 2)      GAZ-08

UniFit Right Angle Connector  
(PKT 3)      GAZ-04-90



### Quartz Adaptor

21-809-3772



### Quartz Adaptor with Aux Port

21-809-3771



## Spray Chambers for HF

### Tracey TFE

50mL Cyclonic      20-809-3717

Suitable for use with HF.  
Requires Transfer Tube 21-809-3726.



### PuraFlex Transfer Tube

21-809-3726



## Torches

### Quartz Torch

SpectroBlue EOP	31-808-3394
SpectroBlue TI	31-808-3468



### Injectors for Quartz Torch

Quartz Injector 1.2mm	31-808-3410
Quartz Injector 1.8mm	31-808-3391
Quartz Injector 2.0mm	31-808-3411
Quartz Injector 3.0mm	31-808-3412
Alumina Injector 1.8mm	31-808-3409
Alumina Injector 2.5mm	31-808-3390



## Common Accessories

### Elegra Argon Humidifier

Elegra Argon Humidifier	70-803-1269
Elegra Stand	70-803-1581




### Guardian In-Line Filters

Guardian DC in-line non-return gas filter	70-803-1942
Guardian in-line sample filter	70-803-1108
Eluo Adaptor for Inline Filter	70-803-1160



### Trident Internal Standard Kits

Trident CT Standard Kit	60-703-1179
Standard Kit for HF	60-808-1150



### Magnifier Inspection Tool 70-803-1923



## Autosampler Probes (See page 22-23 for complete list of probes)

### For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm	
Probe Connecting Line (Red)	70-803-1803
Polyimide sheathed Autosampler Probe	
0.5mm ID	60-808-1186L
PTFE Encapsulated Carbon Fibre Probe	
0.5mm ID	70-803-0784
PTFE Encapsulated Carbon Fibre Probe	
0.75mm ID	70-803-1443
PTFE Encapsulated Carbon Fibre Probe	
1.0mm ID	70-803-0793



### For Cetac ASX-7400, 7600 Series

Guardian Probe 0.75mm	
Probe Connecting Line (Red)	70-803-2940

